


Get to know your

MED-EL RONDO 3 Speech Processor



The processor goes on the charging unit with the flat side facing upwards. The orange light will come on to indicate that the device is charging. It will go off once it is fully charged (after approximately 4 hours).

TOP TIPS

- Your user manual or the MED-EL website can help you to troubleshoot problems: <http://www.medel.com/uk/troubleshooting/>
- Or visit the AIS website for helpful information and troubleshooting advice.
- If you need further help, contact the repairs service.

Get to know your

MED-EL Fine Tuner Echo



TOP TIPS

- If you have a new or replacement processor, you will need to pair it with the fine tuner echo. To do this, place the processor on the back of the fine tuner and switch the processor off and back on again. To pair two processors to one fine tuner, repeat the above step with your second processor.
- The battery in the fine tuner should last for a long time. We do keep spares so please request one if your battery does run out.
- For bilateral users, decide which processor(s) you wish to change and push the appropriate button on the bottom row before pushing a different button(s).